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1. Poynard T, Morra R, Ingiliz P, Imbert-Bismut F, Thabut D, Messous D, Munteanu M, Massard J, Benhamou Y, Ratziu V. Biomarkers of liver fibrosis. *Adv Clin Chem.* 2008;46:131-60. Review.
2. Halfon P, Munteanu M, Poynard T. FibroTest-ActiTest as a non-invasive marker of liver fibrosis. *Gastroenterol Clin Biol.* 2008;32:22-39.
3. Poynard T, Morra R, Halfon P, Castera L, Ratziu V, Imbert-Bismut F, Naveau S, Thabut D, Lebrech D, Zoulim F, Bourliere M, Cacoub P, Messous D, Munteanu M, de Ledinghen V. Meta-analyses of FibroTest diagnostic value in chronic liver disease. *BMC Gastroenterol.* 2007;7:40.
4. Halfon P, Imbert-Bismut F, Messous D, Antoniotti G, Benchetrit D, Cart-Lamy P, Delaporte G, Doutheau D, Klump T, Sala M, Thibaud D, Trepo E, Thabut D, Myers RP, Poynard T. A prospective assessment of the inter-laboratory variability of biochemical markers of fibrosis (FibroTest) and activity (ActiTest) in patients with chronic liver disease. *Comp Hepatol.* 2002;1:3.
5. Imbert-Bismut F, Messous D, Thibaut V, Myers RB, Piton A, Thabut D, Devers L, Hainque B, Mercadier A, Poynard T. Intralaboratory analytical variability of biochemical markers of fibrosis (FibroTest) and activity (ActiTest) and reference ranges in healthy blood donors. *Clin Chem Lab Med.* 2004;42:323-33.
6. Morra R, Munteanu M, Imbert-Bismut F, Messous D, Ratziu V, Poynard T. FibroMax: towards a new universal biomarker of liver disease? *Expert Rev Mol Diagn.* 2007;7:481-90. Review.
7. Munteanu M, Messous D, Thabut D, Imbert-Bismut F, Jouys M, Massard J, Piton A, Bonyhay L, Ratziu V, Hainque B, Poynard T. Intra-individual fasting versus postprandial variation of biochemical markers of liver fibrosis (FibroTest) and activity (ActiTest). *Comp Hepatol.* 2004;3:3.
8. Messous D, Munteanu M, Morra R, Piton A, Poynard T, Hainque B, Imbert-Bismut F. Time impact of serum storage temperature on stability of important liver enzymes and proteins, and on FibroTest-ActiTest results. *Clinica Chimica Acta.* 2005;335:400.
9. Rosenthal-Allieri MA, Peritore M, Tran A, Halfon P, Benzakena S, Bernada A. Analytical variability of the FibroTest proteins. *Clin Biochem* 2005;38:473-8.
10. Rosenthal-Allieri MA, Tran A, Halfon P, Imbert-Bismut F, Munteanu M, Messous D, Peritore ML, Poynard T, Bernard A. Optimal correlation between different instruments for FibroTest-ActiTest protein measurement in patients with chronic hepatitis C. *Gastroenterol Clin Biol.* 2007;31:815-21.
11. Piton A, Messous D, Imbert-Bismut F, Berges J, Munteanu M, Poynard T, Hainque B. Alpha2-macroglobulin immunoturbidimetric assays (DakoCytomation reagents) on Roche Diagnostic analyzers (Modular P, Cobas Integra). Application to FibroTest-ActiTest (FT-AT). *Ann Biol Clin.* 2005;63:385-95.
12. Imbert-Bismut F, Messous D, Raoult A, Poynard T, Bertrand JJ, Marie PA, Louis V, Audy C, Thouy JM, Hainque B, Piton A. Results transferability on RXL, ARX, X-Pand, BN2 (Dade Behring) and modular DP (Roche Diagnostics) analyzers: application to component assays of FibroTest and ActiTest. *Ann Biol Clin.* 2005;63:305-13.
13. Bienvenu J, Doche C, Later R, Manceaux J.C, Pontet E, Pressac M, Vassault A, Dumont G. Améliorations apportées par le matériau de référence CRM 470 dans la standardisation du dosage des protéines sériques. *Ann Bio Clin.* 1997;55:37-40.
14. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J, Blaabjerg O, Blirup- Jensen S, Carlstrom A, Petersen PH, Johnson AM, Milford- Ward A, Ritchie RF, Svendsen PJ, Whicher J. Consensus of a group of professional societies and diagnostic companies for interim reference

- ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP Reference Material (CRM 470). International Federation of Clinical Chemistry. Community Bureau of Reference of the commission of the European Communities. College of American Pathologists. Eur J Clin Chem Clin Biochem. 1996;34:517-20.
15. Albers J.J, Marcovina S.M, and H. Kennedy International Federation of Clinical Chemistry Standardization Project for Measurement of Apolipoproteins A1 and B. II Evaluation and Selection of Candidate Reference Materials. Clin Chem. 1992;38:658-62.
 16. Myara A, Ferard G, Guéchet J, Imbert-Bismut F, Lasnier E, Piton A, Voitot H, Dhumeaux D. Calibration to achieve standardization in enzymology, example of ALT. Clin Chem Lab Med. 2001;39, Special Supplement pp S1-S448.
 17. Groupe de Travail SFBC (Lessinger J.M, Doffoel S, Ferard G, Grafmeyer D, Schiele F, Vassault A, Vialle A) -Assurance de Qualité en Enzymologie Clinique- et Hépatites-(Myara A, Guéchet J, Imbert-Bismut F, Lasnier E, Piton A, Voitot H). Le calibrage en enzymologie clinique : principe, conditions d'application et résultats. Ann Biol Clin. 2002;60:281-6.
 18. Ferard G, Piton A, Messous D, Imbert-Bismut F, Frairi A, Poynard T, Lessinger JM. Intermethod calibration of alanine aminotransferase (ALT) and gamma-glutamyltransferase (GGT) results: application to FibroTest and ActiTest scores. Clin Chem Lab Med. 2006; 44: 400-6.
 19. Schumann G, Bonora R, Ceriotti F, Féraud G, Ferrero CA, Franck PF, Gella FJ, Hoelzel W, Jørgensen PJ, Kanno T, Kessner A, Klauke R, Kristiansen N, Lessinger JM, Linsinger TP, Misaki H, Panteghini M, Pauwels J, Schiele F, Schimmel HG, Weidemann G, Siekmann L; International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of catalytic concentration of gamma-glutamyltransferase. Clin Chem Lab Med. 2002;40:734-8.
 20. Schumann G, Bonora R, Ceriotti F, Féraud G, Ferrero CA, Franck PF, Gella FJ, Hoelzel W, Jørgensen PJ, Kanno T, Kessner A, Klauke R, Kristiansen N, Lessinger JM, Linsinger TP, Misaki H, Panteghini M, Pauwels J, Schiele F, Schimmel HG, Weidemann G, Siekmann L; International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 4. Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. Clin Chem Lab Med. 2002;40:718-724.
 21. Schumann G, Bonora R, Ceriotti F, Féraud G, Ferrero CA, Franck PF, Gella FJ, Hoelzel W, Jørgensen PJ, Kanno T, Kessner A, Klauke R, Kristiansen N, Lessinger JM, Linsinger TP, Misaki H, Panteghini M, Pauwels J, Schiele F, Schimmel HG, Weidemann G, Siekmann L; International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 5. Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med. 2002;40:725-731.
 22. Féraud G, Imbert-Bismut F, Messous D, Piton A, Abella A, Burnat P, Hainque B, Glasser N, Lessinger JM. Influence of pyridoxal phosphate in measuring aminotransferases activities in patients with viral hepatitis. Ann Biol Clin. 2004;62:717-20. French.
 23. Jendrassik Grof. Vereinfachte Photometrische Methoden zur Bestimmung des Blutbilirubins. Biochem Z. 1938;297:8189.
 24. Malloy HT, Evelyn KA. The determination of bilirubin with photoelectric colorimeter J Biol Chem. 1937;19:984-93.

25. Wahlefeld AW, Herz G, Bernt E. Modification of the Malloy-Evelyn method for a simple, reliable determination of total bilirubin in serum. *Scand J Clin Lab Invest.* 1972;29,supp 26, Abstract 11-12.
26. Doumas BT, Perry BW, Sasse EA, Straumfjord JV Jr. Standardization in bilirubin assays: evaluation of selected methods and stability of bilirubin solutions. *Clin Chem.* 1973;19:984-93.
27. Trinder P. *Ann Clin Biochem.* 1969;6:24.
28. Peterson JL and Young DS. *Anal Biochemistry.* 1958;23:301.
29. Thierry Poynard, Mona Munteanu, Olivier Deckmyn, Yen Ngo, Fabienne Drane, Djamila Messous, Jean Marie Castille, Chantal Housset, Vlad Ratziu, Françoise Imbert-Bismut. Applicability and precautions of use of liver injury biomarker FibroTest. A reappraisal at 7 years of age. *BMC Gastroenterol.* 2011;11:39.
30. Bedossa P, Dargère D, Paradis V. Sampling variability of liver fibrosis in chronic hepatitis C. *Hepatology.* 2003;38:1449-57.
31. Poynard T, Munteanu M, Imbert-Bismut F, Charlotte F, Thabut D, Le Calvez S, Messous D, Thibault V, Benhamou Y, Moussalli J, Ratziu V. Prospective analysis of discordant results between biochemical markers and biopsy in patients with chronic hepatitis C. *Clin Chem.* 2004;50:1344-55.
32. Costa JM, Telehin D, Munteanu M, Kobryn T, Ngo Y, Thibault V, Joseph M, Ratziu V, Benhamou Y, Koz'ko V, Dubins'ka G, Poveda JD, Poynard T. HCV-GenoFibrotest: a combination of viral, liver and genomic (IL28b, ITPA, UGT1A1) biomarkers for predicting treatment response in patients with chronic hepatitis C. *Clin Res Hepatol Gastroenterol.* 2011;35:204-13.
33. Poynard T, de Ledinghen V, Zarski JP, Stanciu C, Munteanu M, Vergniol J, France J, Trifan A, Lenaour G, Vaillant JC, Ratziu V, Charlotte F; Fibrosis-TAGS group. Performances of Elasto-FibroTest(®), a combination between FibroTest(®) and liver stiffness measurements for assessing the stage of liver fibrosis in patients with chronic hepatitis C. *Clin Res Hepatol Gastroenterol.* 2012;36:455-63.
34. Deckmyn O. BioPredictive 2010. SteatoScreen internal evaluation.
35. Poynard T, Peta V, Munteanu M, Charlotte F, Ngo Y, Ngo A, Perazzo H, Deckmyn O, Pais R, Mathurin P, Myers R, Loomba R, Ratziu V; FLIP consortium, the FibroFrance-CPAM group, the FibroFrance-Obese group, and the Selonsertib group. The diagnostic performance of a simplified blood test (SteatoTest-2) for the prediction of liver steatosis. *Eur J Gastroenterol Hepatol.* 2019;31:393-402.
36. Poynard T, Munteanu M, Charlotte F, Perazzo H, Ngo Y, Deckmyn O, Pais R, Merrouche W, de Ledinghen V, Mathurin P, Ratziu V; FLIP Consortium, the FibroFrance-CPAM Group, the FibroFrance-Obese Group. Diagnostic performance of a new noninvasive test for nonalcoholic steatohepatitis using a simplified histological reference. *Eur J Gastroenterol Hepatol.* 2018;30:569-77.
37. Gossens K, Van Uytvanghe K, Thienpont LM. Trueness and comparability assessment of widely used assays for 5 common enzymes and 3 electrolytes. *Clin Chim Acta.* 2015;442:44-5.
38. Infusino I, Frusciante E, Braga F, Panteghini M. Progress and impact of enzyme measurement standardization. *Clin Chem Lab Med.* 2017;55:334-40.
39. Varaut A, Fontaine H, Serpaggi J, Verkarre V, Vallet-Pichard A, Nalpas B, Imbertbismuth F, Lebray P, Pol S. Diagnostic accuracy of the fibrotest in hemodialysis and renal transplant patients with chronic hepatitis C virus. *Transplantation.* 2005;80:1550-5.
40. Alric L, Kamar N, Bonnet D, Danjoux M, Abravanel F, Lauwers-Cances V, Rostaing L. Comparison of liver stiffness, fibrotest and liver biopsy for assessment of liver fibrosis in kidney-transplant patients with chronic viral hepatitis. *Transpl Int.* 2009;22:568-73.
41. Ratziu V, Charlotte F, Heurtier A, Gombert S, Giral P, Bruckert E, Grimaldi A, Capron F, Poynard T; LIDO Study Group. Sampling variability of liver biopsy in nonalcoholic fatty liver disease. *Gastroenterology.* 2005;128:1898-906.

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